

How to Do a SARS-CoV-2 Antigen Rapid Diagnostic Test

Sample collection, testing and interpretation of SARS-CoV-2 Antigen RDTs

READ THE MANUFACTURERS' INSTRUCTIONS CAREFULLY BEFORE YOU BEGIN


Preparation for sample collection

- A. New (unopened) individually wrapped sterile nasopharyngeal swab
- B. Personal protective equipment-Gloves, gown, eye protection or face shield & N95 respirator
- C. Pen for marking or labelling






Preparation for SARS-CoV-2 Antigen testing

- A. New (unopened) SARS-CoV-2 Antigen Test Device
- B. Timer
- C. Extraction buffer tube and nozzle
- D. Medical waste container
- E. Disinfectant (bleach and alcohol)






Sample collection

- 1** Insert a sterile nasopharyngeal swab into the nasal cavity of the patient, reaching the surface of the posterior nasopharynx
- 2** Swab over the surface of the posterior nasopharynx, rotating the swab 3-4 times to ensure a good sample. Leaving the swab in the nasal cavity for a few seconds will ensure absorption of the nasal secretions
- 3** Withdraw the sterile swab from the nasal cavity



SARS-CoV-2 Antigen testing

- 4** Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab
- 5** Remove the swab while squeezing the sides of the tube to extract the liquid from the swab
- 6** Press the nozzle cap tightly onto the tube
- 7** Apply drops of extracted specimen to the specimen well of the test device. Add the exact number of drops specified by the manufacturer
- 8** Read and record the test result after the specified period. Report the test result
- 9** Dispose of the gloves, test device, swab and extraction buffer tube in the medical waste container. Disinfect the workstation with bleach and alcohol



Interpretation of the test result

Positive

A line in "C" AND a line in "T" means SARS-CoV-2 is DETECTED



Negative

A line in "C" and NO LINE in "T" means SARS-CoV-2 is NOT DETECTED



Invalid

NO LINE in "C" and a line or no line in "T" means the test is INVALID. Repeat the test using a new (unopened) SARS-CoV-2 Antigen test device and a new sample

